Interpreting results

Data generated using Advanced Instruments standards and reference solutions may be analyzed according to the accuracy and precision specifications of the instrument (see *Performance characteristics - Osmo1 Single-sample Micro-osmometer User's Guide*).

Laboratories may choose to employ one, two, or three standard deviations (SD) for accuracy based on what is relevant in their laboratories. For normally distributed data, approximately 68% of the individual data values will fall within one standard deviation of the mean, approximately 95% within two standard deviations, and approximately 99.7% within three standard deviations.

Advanced Instruments Standards and Reference Solutions	Precision	Accuracy (mOsm/kg H ₂ O)		
		1 SD	2 SD	3 SD
Clinitrol™ 290 Reference Solution	Standard deviation ≤2 mOsm/kg H ₂ O	288-292	286-294	284-296
50 mOsm/kg Calibration Standard	Standard deviation ≤2 mOsm/kg H ₂ O	48-52	46-54	44-56
100 mOsm/kg Calibration Standard	Standard deviation ≤2 mOsm/kg H ₂ O	98-102	96-104	94-106
200 mOsm/kg Calibration Standard	Standard deviation ≤2 mOsm/kg H ₂ O	198-202	196-204	194-206
400 mOsm/kg Calibration Standard	Standard deviation ≤2 mOsm/kg H ₂ O	398-402	396-404	394-406
500 mOsm/kg Calibration Standard	Coefficient of variation ≤0.5%	497-503	494-506	491-509
850 mOsm/kg Calibration Standard	Coefficient of variation ≤0.5%	845-855	840-860	835-865
900 mOsm/kg Calibration Standard	Coefficient of variation ≤0.5%	895-905	890-910	885-915
1000 mOsm/kg Calibration Standard	Coefficient of variation ≤0.5%	995-1005	990-1010	985-1015
1500 mOsm/kg Calibration Standard	Coefficient of variation ≤1%	1485-1515	1470-1530	1455-1545
2000 mOsm/kg Calibration Standard	Coefficient of variation ≤1%	1980-2020	1960-2040	1940-2060

Advanced Instruments Controls	Expected Levels	Expected Ranges
Protinol™ Protein-Based Controls	240 mOsm/kg H ₂ O 280 mOsm/kg H ₂ O 320 mOsm/kg H ₂ O	233-247 mOsm/kg H ₂ O 273-287 mOsm/kg H ₂ O 313-327 mOsm/kg H ₂ O
Renol™ Urine Osmolality Controls	300 mOsm/kg H ₂ O 800 mOsm/kg H ₂ O	290-310 mOsm/kg H ₂ O 790-810 mOsm/kg H ₂ O

Protinol and Renol specifications are presented as a range. This range was established based of stability data and approximates a three standard deviation range of results.

Two Technology Way | Norwood, MA 02062

800-225-4034 | +1-781-320-9000

aicompanies.com | info@aicompanies.com

Technical Service Advanced Instruments and worldwide distributor network provides 24/7 comprehensive customer service and technical support.

This document is copyrighted by Advanced Instruments with all rights reserved. Under copyright laws, this guide may not be reproduced in any form, in whole or part, without the prior written consent of Advanced Instruments.

© 2020 Advanced Instruments. Osmo1 is a trademark of Advanced Instruments. All other trademarks are the property of their respective companies.

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