## Total nucleated cell (TNC) reagent

**REF** 111012 Three (3) 1.8 mL vials

### Intended use

For use on the Advanced Instruments GloCyte Automated Cell Counter System only. The Total Nucleated Cell (TNC) Reagent stains total nucleated cells for quantitative enumeration of TNCs in human cerebrospinal fluid (CSF) specimens.

### Principles of the method

A detergent in the reagent permeabilizes the cell membrane of nucleated cells, and a fluorescent dye binds to the nucleic acid present in the cells, fluorescently labeling the nucleated cells. A 30 µL sample of neat or diluted CSF is stained with 30 µL of TNC Reagent in a tube. A 30 µL or 10 µL sample of the stained CSF specimen is dispensed onto a test cartridge, where the cells are captured on the membrane using the vacuum station. The test cartridge is placed into the GloCyte instrument, and the stained TNCs are imaged and enumerated using fluorescence microscopy and digital image analysis. Refer to the *User's Guide*<sup>1</sup> for detailed instructions.

### Product composition

The TNC Reagent contains < 3 µg/mL propidium iodide in a buffered solution with < 0.1% w/v sodium azide as preservative.

### Preparation

The TNC Reagent is a liquid, ready-to-use.

### Warnings and precautions

TNC Reagent should be handled according to established Good Laboratory Practices and Universal Precautions. Dispose of materials according to your institution's practices in a safe and acceptable manner that is in compliance with all local, regional, national and/or international regulations.

**Caution:** Avoid contact with metals because sodium azide may react with lead or copper plumbing to form highly explosive metal azides; flush with copious amounts of water when pouring dilute solutions down the drain to prevent buildup. Propidium iodide is a mutagen; handle with care.

### Storage, stability and handling

	Storage	Stability
Unopened	2-8 °C (36-46 °F)	Refer to vial label
Opened	2-25 °C (36-77 °F)	2 weeks

- Tightly seal the vials immediately after use.
- Do not mix vigorously.
- Do not wipe excess liquid from external vial threads.
- Always store in original vials.
- Do not freeze.

### Specimen collection and handling

Collect CSF in a sterile CSF collection tube. The CSF specimen should be at room temperature, 15-25 °C (59-77 °F). Testing should be completed as soon as possible, as cellular degeneration in CSF can begin within one hour of collection<sup>2</sup>. Hand-mix the specimen just prior to testing by inverting the specimen tube 10-15 times.

### Instructions for use

Please refer to the *User's Guide*<sup>1</sup>.

### Interfering substances





- Cellular debris from hemolysis may cause erroneous results for the GloCyte TNC assay.
- Extremely high levels of some bacteria (e.g., *Neisseria lactamica*) at >10<sup>7</sup> CFU/mL may cause erroneous results for the GloCyte TNC assay.

Nucleated red blood cells are counted by the GloCyte TNC assay. Erroneous results can occur from adverse storage conditions or use of expired materials.

### For sales and service

Contact your Advanced Instruments Distributor.

### Prescription use only

Two Technology Way / 781-320-9000  
Norwood, Massachusetts 02062, USA  
800-225-4034 Fax: 781-320-8181  
aicompanies.com

© 2021 Advanced Instruments

GloCyte is a registered trademark of Advanced Instruments

<sup>1</sup> GloCyte Automated Cell Counter System User's Guide (p/n 111005UG)  
<sup>2</sup> CLSI H56-A, Body Fluid Analysis for Cellular Composition; Approved Guideline, 2006.

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WMDE B.V.