Red blood cell (RBC) reagent

REF 111011 Three (3) 1.8 mL vials

Intended use

For use on the Advanced Instruments GloCyte Automated Cell Counter System only. The Red Blood Cell (RBC) Reagent stains erythrocytes for quantitative enumeration of RBCs in human cerebrospinal fluid (CSF) specimens.

Principles of the method

An anti-human RBC antibody conjugated to a fluorescent dye in the reagent binds to the membranes of erythrocytes, fluorescently labeling the RBCs. A 30 µL sample of neat or diluted CSF is stained with 30 µL of RBC Reagent in a tube for 3 minutes. After incubation, a 30 µL or 10 µL sample of the stained specimen is dispensed onto a test cartridge, where the cells are captured on the membrane using the vacuum station. The test cartridge is placed into the GloCyte instrument, and the stained RBCs are imaged and enumerated using fluorescent microscopy and digital image analysis. Refer to the *User's Guide*¹ for detailed instructions.

Product composition

The RBC Reagent contains < 2 µg/mL mouse anti-human RBC antibody with fluorochrome in a buffered solution with < 1% w/v bovine serum albumin and < 0.1% w/v sodium azide as preservative.

Preparation

The RBC Reagent is a liquid, ready-to-use.

Warnings and precautions

RBC Reagent should be handled according to established Good Laboratory Practices and Universal Precautions. Dispose of materials according to your institution's practices in a safe and acceptable manner that is in compliance with all local, regional, national, and/or international regulations.

Caution: Avoid contact with metals; sodium azide may react with lead or copper plumbing to form highly explosive metal azides; flush with copious amounts of water when pouring dilute solutions down the drain to prevent buildup.

Storage, stability, and handling

	Storage	Stability
Unopened	2 - 8 °C (36 - 46 °F)	Refer to vial label
Opened	2 - 25 °C (36 - 77 °F)	2 weeks

- Tightly seal the vials immediately after use.
- Do not mix vigorously.
- Do not wipe excess liquid from external vial threads.
- Always store in original vials.
- Do not freeze.

Specimen collection and handling

Collect CSF in a sterile CSF collection tube. The CSF specimen should be stored at room temperature 15 - 25 °C (59 - 77 °F). Testing should be completed as soon as possible, as cellular degeneration in CSF can begin within one hour of collection². Hand-mix the specimen just prior to testing by inverting the specimen tube 10-15 times.

Instructions for use

Please refer to the *User's Guide*¹.

Interfering substances





RBC fragments and cellular debris from hemolysis may cause erroneous results for the GloCyte RBC assay.

Erroneous results can occur from adverse storage conditions or use of expired materials.

For sales and service

Contact your Advanced Instruments Distributor.

Prescription use only

Two Technology Way / 781-320-9000
Norwood, Massachusetts 02062, USA
800-225-4034 Fax: 781-320-8181
aicompanies.com

© 2021 Advanced Instruments

GloCyte is a registered trademark of Advanced Instruments

¹ GloCyte Automated Cell Counter System User's Guide (p/n 111005UG)
² CLSI H56-A. Body Fluid Analysis for Cellular Composition; Approved Guideline. 2006.

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