Interpreting results

Data generated using Advanced Instruments standards and reference solutions may be analyzed according to the accuracy and precision specifications of the instrument (see *Performance characteristics* in the *OsmoPRO Multi-sample Micro-osmometer User's Guide* (p/n 112005UG)).

Laboratories may choose to employ one, two, or three standard deviations (SD) for accuracy based on what is relevant in their laboratories. For normally distributed data, approximately 68% of the individual data values will fall within one standard deviation of the mean, approximately 95% within two standard deviations, and approximately 99.7% within three standard deviations.

Advanced Instruments Standards and Reference Solutions	Precision	Accuracy (mOsm/kg H ₂ O)		
		1 SD	2 SD	3 SD
Clintrol™ 290 Reference Solution	Standard deviation ≤ 3 mOsm/kg H ₂ O	287-293	284-296	281-299
50 mOsm/kg Calibration Standard	Standard deviation ≤ 3 mOsm/kg H ₂ O	47-53	44-56	41-59
100 mOsm/kg Calibration Standard	Standard deviation ≤ 3 mOsm/kg H ₂ O	97-103	94-106	91-109
200 mOsm/kg Calibration Standard	Standard deviation ≤ 3 mOsm/kg H ₂ O	197-203	194-206	191-209
400 mOsm/kg Calibration Standard	Standard deviation ≤ 3 mOsm/kg H ₂ O	397-403	394-406	391-409
500 mOsm/kg Calibration Standard	Coefficient of variation ≤ 0.75%	496-504	492-508	488-512
850 mOsm/kg Calibration Standard	Coefficient of variation ≤ 0.75%	843 -857	836-864	829-871
900 mOsm/kg Calibration Standard	Coefficient of variation ≤ 0.75%	893-907	886-914	879-921
1000 mOsm/kg Calibration Standard	Coefficient of variation ≤ 0.75%	992-1008	984-1016	976-1024
1500 mOsm/kg Calibration Standard	Coefficient of variation ≤ 0.75%	1488-1512	1476-1524	1464-1536
2000 mOsm/kg Calibration Standard	Coefficient of variation ≤ 0.75%	1985-2015	1970-2030	1955-2045

Advanced Instruments Controls	Expected Levels	Expected Ranges
Protinol™ Protein-Based Controls	240 mOsm/kg H ₂ O 280 mOsm/kg H ₂ O 320 mOsm/kg H ₂ O	233-247 mOsm/kg H ₂ O 273-287 mOsm/kg H ₂ O 313-327 mOsm/kg H ₂ O
Renol™ Urine Osmolality Controls	300 mOsm/kg H ₂ O 800 mOsm/kg H ₂ O	290-310 mOsm/kg H ₂ O 790-810 mOsm/kg H ₂ O

Protinol and Renol specifications are presented as a range. This range was established based on stability data and approximates a three standard deviation range of results.

Two Technology Way / 781-320-9000
Norwood, Massachusetts 02062, USA
800-225-4034 Fax: 781-320-8181
aicompanies.com

Advanced Instruments and worldwide distributor network provides comprehensive customer service and technical support.

This document is copyrighted by Advanced Instruments with all rights reserved. Under copyright laws, this guide may not be reproduced in any form, in whole or part, without the prior written consent of Advanced Instruments.

© 2021 Advanced Instruments. OsmoPRO is a trademark of Advanced Instruments. All other trademarks are the property of their respective companies.

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