



Certificate of Analysis

Product: Broad Spectrum Microtab II – Bulk


Lot No.: 10812

Item Code: 954

Date of Manufacture: August 2012

<u>TEST</u>	<u>SPECIFICATION</u>	<u>RESULT</u>
Description	A small, round, mottled pink, mauve biconvex tablet	Complies
Average Weight	About 18 mg	18.02
Diameter	About 3 mm	3.01
Solution Time (pasteurized milk)	Not more than 30 minutes	14
Solution Appearance (milk)	No specks visible	Complies
Bronopol Content	7.0 to 9.0 mg/tablet	7.5
Natamycin Content	0.25 to 0.35 mg/tablet	.26

Results submitted by an outside laboratory.



Robert Soares
Director of Quality
D&F Control Systems, Inc.
September 20, 2012