



Certificate of Analysis

Product: Broad Spectrum Microtab II – Bulk

Lot No.: K40713

Item Code: 954

Date of Manufacture: August 2013

Expiration Date: August 2017

<u>TEST</u>	<u>SPECIFICATION</u>	<u>RESULT</u>
Description	A small, round, mottled pink, mauve biconvex tablet	Complies
Average Weight	About 18 mg	18.37mg
Diameter	About 3 mm	3.00mm
Solution Time (pasteurized milk)	Not more than 30 minutes	13 min
Solution Appearance (milk)	No specks visible	Complies
Bronopol Content	7.0 to 9.0 mg/tablet	7.60
Natamycin Content	0.25 to 0.35 mg/tablet	0.27

Results submitted by an outside laboratory.



Meredith Pesta
Laboratory Services Manager
D&F Control Systems, Inc.
August 19, 2013